



Sample Preservation Chart

Fresno Lab
687 N. Laverne Ave..
Fresno, CA 93727
P 559.497.2888

Sacramento Lab
9332 Tech Center Dr. Ste. 100
Sacramento CA 95826
P 916.853.9293

San Bernardino Lab
350 E. Commercial Rd. #110
San Bernardino, CA 92408
P 909.796.2059

Vancouver Lab
2517 East Evergreen Blvd.
Vancouver, WA 98661
P 360.750.0055

Additional Drop off Locations in California: Visalia, Bakersfield, Livermore & San Jose Analytical Services provided in CA, NV, OR, WA, & HI

www.bskassociates.com

Organic Analyses

	Method	Analyte Groups	Hold Time	Container	Preservation Requirements	
Drinking Water	EPA 524.2	Trihalomethanes	14 days	3 x 40mL AG VOAs	Na ₂ S ₂ O ₃ + TB	≤ 6°C
		VOCs—Raw Water	14 days	3 x 40mL VOAs	HCl + TB	pH ≤ 2 ≤ 6°C
	EPA 524.2 / CA DHS	VOCs—Finished Water & 1,2,3-Trichloropropane	14 days	250 mL AG & x 40mL VOAs	3 Ascorbic (AG) & HCl (VOAs) + TB	pH ≤ 2 ≤ 6°C
	524-TCP	TCP	14 days	3 x 40mL AG VOAs	HCL + TB	pH ≤ 2 ≤ 6°C
	EPA 504	EDB, DBCP	14 days	3 x 40mL VOAs	Na ₂ S ₂ O ₃ + TB	≤ 6°C
	EPA 505	Organochlorine Pesticides	7 days	3 x 40mL VOAs	Na ₂ S ₂ O ₃	≤ 6°C
	EPA 515	Chlorinated Herbiidices	14 days	250mL AG	Na ₂ SO ₃	≤ 6°C
	EPA 525	Semivolatile Organics	14 days	2 x 1L AG	Ascorbic, EDTA, KH ₂ Ct	≤ 6°C
	EPA 531	Carbamates	28 days	1 x 40mL VOA	MCAA Buffer + Na ₂ S ₂ O ₃	pH = 3 ≤ 6°C
	EPA 537.1	PFAS	14 days	3 X 250ml (P) (FB Required)	Trizma	pH 6-8 ≤ 6°C
EPA 533	PFAS	28 days	3 X 250ml (P) (FB Required)	Ammonium Acetate	pH 6-8 ≤ 6°C	
EPA 547	Glyphosate	14 days	1 x 40mL VOA	Na ₂ S ₂ O ₃	≤ 6°C	
EPA 548	Endothall	7 days	250mL AG	Na ₂ S ₂ O ₃	≤ 6°C	
EPA 549	Diquat/Paraquat	7 days	1L Brown Plastic	Na ₂ S ₂ O ₃	≤ 6°C	
EPA 552	Haloacetic Acids	14 days	250mL AG	NH ₄ Cl	≤ 6°C	
EPA 608	OCI Pesticides, PCBs	7 days	2 x 1L AG	No Preservative	≤ 6°C	
Wastewater	EPA 624	Volatile Organics	14 days	3 x 40mL VOAs	HCl	pH ≤ 2 ≤ 6°C
		Acrolein, Acrylonitrile	14 days	1 x 40mL VOA	pH 4-5 Buffer	pH 4-5 ≤ 6°C
		2-Chloroethyl vinyl ether	14 days	2 x 40mL VOAs	No Preservative	≤ 6°C
	EPA 625	Semivolatile Organics	7 days	2 x 1L AG	No Preservative	≤ 6°C
	EPA 8321	Carbamates	7 days	1 x 40mL AG VOA	No Preservative	≤ 6°C
RCRA	EPA 8081/8082	OCI Pesticides & PCBs	L=7; S=14 days	2 x 1L AG/ ST or Jar	No Preservative	≤ 6°C
	EPA 8151	Phenoxy Acid Herbicides	L=7; S=14 days	2 x 1L AG/ ST or Jar	No Preservative	≤ 6°C
	EPA 8260	Volatile Organics	14 days	3 x 40mL VOAs	HCl	pH ≤ 2 ≤ 6°C
		Acrolein, Acrylonitrile	7 days	1 x 40mL VOA	pH 4-5 Buffer	pH 4-5 ≤ 6°C
		2-Chloroethyl vinyl ether	7 days	1 x 40mL VOA	No Preservative	≤ 6°C
Hydro	EPA 8270	Semivolatile Organics	L=7; S=14 days	2 x 1L AG/ ST or Jar	No Preservative	≤ 6°C
	GC/MS	BTEX / TPH-G	14 days	3 x 40mL VOAs	HCl	pH ≤ 2 ≤ 6°C
	EPA 1664	Oil & Grease HEM	28 days	2 x 1L AG	HCl	pH ≤ 2 ≤ 6°C
	EPA 8015	Diesel	L=7; S=14 days	2 x 1L AG / ST or Jar	HCl	pH ≤ 2 ≤ 6°C

ST=Soil Tube AG=Amber Glass 1L=1 Liter L=Liquid S=Solid TB=Trip / Travel Blank

	Analyses	Hold Time	Container	Preservation Requirements			
Wet Chemistry	Alkalinity, Bicarbonate, Carbonate	14 days	250mL Plastic	No Preservative	≤ 6°C		
	BOD, CBOD (call for instructions)	48 hours	1L Plastic	No Preservative	≤ 6°C		
	Bromate, Chlorate	28 days	250mL AG	EDA	≤ 6°C		
	Chlorite	14 days	250mL AG	EDA	≤ 6°C		
	Bromide	28 days	250mL Plastic or AG	EDA or No Preserv.			
	Carbon, Total Org. (TOC)	28 days	3 x 40mL VOAs	H ₃ PO ₄	pH ≤ 2	≤ 6°C	
	Carbon, Dissolved Org. (DOC)	28 days	3 x 40mL VOAs	No Preservative		≤ 6°C	
	COD	28 days	250mL AG	H ₂ SO ₄	pH ≤ 2	≤ 6°C	
	Nitrate, Nitrite, Orthophosphate	48 hours	250mL Plastic	No Preservative		≤ 6°C	
	Chlorine, Residual and/or pH	15 minutes	250mL Plastic	No Preservative		≤ 6°C	
	Chromium, Hexavalent (Cr+6) DW	14 days	250mL Plastic	NH ₄ OH/(NH ₄) ₂ SO ₄	pH > 8	≤ 6°C	
	Chromium, Hexavalent (Cr+6) WW	28 days	250mL Plastic	NH ₄ OH/(NH ₄) ₂ SO ₄	pH = 9.3-9.7	≤ 6°C	
	Chromium, Hexavalent (Cr+6) WW	24 hours	250mL Plastic	NH ₄ OH/(NH ₄) ₂ SO ₄	pH = 9-9.5	≤ 6°C	
	Color, Turbidity	48 hours	500mL AG	No Preservative		≤ 6°C	
	Conductivity (EC), Sulfate, Chloride, Fluoride	28 days	250mL Plastic	No Preservative		≤ 6°C	
	Cyanide, Total (CN)	14 days	250mL Plastic	NaOH	pH ≥ 12	≤ 6°C	
	Dissolved Oxygen (DO)	15 minutes	300mL Plastic	Special (ABC) vials			
	TKN, Total Phosphorus, Ammonia (NH ₃)	28 days	250mL Plastic	H ₂ SO ₄	pH ≤ 2	≤ 6°C	
	Odor	24 hours	500mL AG	No Preservative		≤ 6°C	
	Perchlorate	28 days	250mL Plastic	No Preservative			
	Solids (TDS,TSS,TS, Volatile and Fixed)	7 days	1L Plastic	No Preservative		≤ 6°C	
	Settable Matter (SM)	48 hours	1L Plastic	No Preservative		≤ 6°C	
	Surfactants (MBAS)	48 hours	1L Plastic	No Preservative		≤ 6°C	
	Gross Alpha	6 months	2 x 1L Plastics	HNO ₃	pH ≤ 2		
	Metals	Al, Sb, As, Ba, Be, B, Cd, Ca, Cr, Co, Cu, Fe, Pb, Mg, Hg, Mn, Mo, Ni, K, Se, Ag, Na, Sr, Ti, Tl, Sn, U, V, Zn	6 months (Hg 28 days)	250mL or 500mL Plastic	HNO ₃	pH ≤ 2	
		Mercury (Hg) by EPA 245.7 (WW only)	28 days	250mL Polyethylene	HCl	pH ≤ 2	
		Silica	28 days	500mL Plastic	No Preservative		≤ 6°C
		Lead/Copper Rule	14 days	1L Plastic	No Preservative		
	Externals	Asbestos	48 hours	L= 1L Plastic	No Preservative		≤ 6°C
		Dioxin-1613	1 year	2 x 1L AG	No Preservative		≤ 6°C
Gross Alpha or Radium 226/228		6 months	2 x 1L Plastics	HNO ₃	pH ≤ 2		
Radon		48 hours	1 - 40ml VOA	No Preservative			
Phenolics-Low Level		28 days	2 x 1L AG	H ₂ SO ₄	pH ≤ 2	≤ 6°C	
Hydrogen Sulfide-Unionized		48 hours	500mL Preserve, 250ml Plastic- None	ZnAc ₂ + NaOH	pH ≥ 9	≤ 6°C	
Total Sulfide		7 Days	500ml Plastic	ZnAc ₂ + NaOH	pH ≥ 12	≤ 6°C	
Dissolved Sulfide		7 Days	500ml Plastic	NaOH	pH ≥ 12	≤ 6°C	
Sulfite (SO ₃)		15 minutes	250mL Plastic	No Preservative		≤ 6°C	
TPH- Diesel		L=7; S=14 days	3- 40ml VOA ST or jar	No Preservative		≤ 6°C	